Your Best Partner for Lifelong Health
你一生健康的最佳伴侶
Rakan Kongsi Kesihatan Terbaik Sepanjang Hayat Anda

Your Best Partner for Lifelong Health
你一生健康的最佳伴侶
Rakan Kongsi Kesihatan Terbaik Sepanjang Hayat Anda

Vertical Vibration Therapy 垂直律動運動 Terapi Getaran Tegak

Over the years, numerous medical experiments have proved that Whole Body Vertical Vibration Therapy achieve overall health benefits through gravitational effects.

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Whole Body Vertical Vibration Training

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The Evolution of Whole Body Vertical Vibration

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The earliest concept of Whole Body Vibration training began in the year 1900. A Whole Body Vibration chair was invented by Dr John Harvey Kellogg (1852 - 1943) in the United States of America for his patients' treatment. However, the modern Whole Body Vertical Vibration device was invented by the Soviet Union in the 1960s and was used for their astronauts' training to improve muscle strength, bone density, endocrine disorder and etc. It is to resolve the health issues faced by astronauts training under a zero gravity environment for a prolonged period of time and to prolong their time spent living in space.

Under zero gravity condition, an astronaut's muscle will begin to shrink and bone density of spine, hipbone and femur will deteriorate rapidly at an average rate of 1 - 1.6% per month. If left unresolved, an astronaut can spend living only around 300 days in space. However, when the Soviet Union's astronaut created a record by staying in space for 437 days through Vertical Vibration Training, which is far longer than any American astronauts, it received worldwide attention.

In the 1980s, the application of Whole Body Vertical Vibration had expanded to athletes' training. Athletes who had undergone this training method showed significant improvement in their performances, thus it began to be widely researched and applied. After the 21st century, this invention has further expanded its application to other medical fields such as medical rehabilitation, sports medicine and etc.

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
Bones
骨頭
Tulang-tulang

Muscles
肌肉
Otot-otot

Principles of Whole Body Vertical Vibration

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Prinsip Getaran Tegak Seluruh Badan

The concept of Whole Body Vertical Vibration originated from the gravity's influence on both human growth and movement. Gravity works on human vertically from bottom to top, and from the external to the internal to achieve the effects of a whole body exercise. In simpler terms, it changes the body's muscle, bone and nervous system through the mechanism of gravity, which further influences the endocrine and other physiological systems in order to produce various beneficial health effects. It is the most natural, safest, most effective and fastest way to exercise.

The Vertical Vibration stimulant will drive an up and down body motion with the speed and acceleration of sine wave through the vibration platform while not burdening the heart and lungs. It is an interval exercise that delivers vibration waves in constant speed from head to toe, mildly stimulating every part of your body, from cells, blood vessels, muscles, nervous system, endocrine system, bones to the organ systems of the body, achieving the effect of repetitive impact training.

The best exercise method in sports medicine is to maximise relaxation of the body, which induces contraction and stretching of the muscle, thus promoting the body's metabolism. This is because the more relaxed the body is the more the blood vessel dilates; blood flow resistance tends to be the lowest and peripheral circulation tends to be the best, minimising cardiopulmonary load. This exercise method is the safest and most effective, and maximises the exercise's benefits on promoting overall blood flow.

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Safety Instructions of Whole Body Vertical Vibration Training

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Arahan Keselamatan Latihan Getaran Tegak Seluruh Badan

PROHIBITED USERS

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Pregnant women and women who are trying to conceive and a month after delivery, patients who have not recovered from hernia, fractures or wounds, individuals who require worry-free recuperation, and children below the age of 6.

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Individuals who are suited for LOW-SPEED VIBRATION

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Individu yang bersesuaian dengan getaran berkelajuan rendah

Women during menstruation, those who can undergo rehabilitation and walking exercise as advised by medical professionals, and individuals with the following physical conditions: those with artificial organs in the body, severe brain and cardiovascular diseases, having recovered from spinal injury, stroke, arthritis, aged and weak. However, for those who have just started, please start with a low speed setting that one is most comfortable with.

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Individual with SPINAL INJURY, OSTEOPOROSIS

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Do not stand straight or sit directly onto the vibration device. Stand with body slightly tilted to the front and hold the table or chair with the hands, and start with low speed or at one's most comfortable speed level. Alternatively, adopt a slight squat position, or sit on a chair with feet on the vibration device to carry out lower body vibration training.

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Individual with HERNIATED DISC

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Individu dengan cakera hernia

Do not stand straight or sit directly on the vibration device. Stand on the vibration device with hands holding onto a table or chair, and lean the body towards the table or chair, at the most comfortable speed. This will not only avoid pressure on the nerves imposed by the bulging disc, but also relaxes and alleviates discomfort on the lower back too. Those who have recovered can adopt a slight squat position and carry out training at a low speed.

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Whole Body Vertical Vibration Training Guide

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Kaedah Kecergasan Getaran Tegak Seluruh Badan

Before Training 律動前 Sebelum Getaran

Drink 200-300ml warm water. This will dilute the blood's concentration; reduce blood flow resistance and blood pressure.

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
30 minutes each time
At least twice a day,
morning and evening

After Training 律動後 Selepas Getaran

Drink 200-300ml warm water. This will promote metabolism and replenish the loss of fluids during vibration.

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Timing of training: Use 1 hour before or after meals

Speed of vibration: For best effects, use at individual's most comfortable speed

Duration of training: 15 minutes each time for the first week, 30 minutes each time from the second week onwards. Best exercise frequency is at least 2 times daily, morning and evening.

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Use 1 hour after meals, do not eat or drink anything other than warm water within 2 hours after training. For a better skin firming effect, drink non-GMO sugar-free warm soy milk. For those with low blood sugar, a glass of warm soy milk will help if feeling hungry after training.

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Have a meal following a 5 minutes rest after training will help you become stronger physically month after month.

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Possible Reactions during Whole Body Vertical Vibration

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Tindak Balas yang Dijanjikan semasa Getaran Tegak Seluruh Badan

ITCHINESS 癢 KEGATALAN

Itching sensation at some parts of the body during vibration training is due to the acceleration of circulation in the body. This itching sensation which is from the inner body will cease with continuous use.

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NUMBNESS 麻 KEKEBASAN

Blood flow will accelerate during vibration training. The continuous impact of the blood when encountering an obstructed area will bring a sense of numbness, which is a good sign. If the obstruction is not cleared, congestion will occur over time causing pain. With continuous vibration training, the condition will improve.

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SORENESS 酸 KERENGAAN

Joints, waist, neck, back, hands or muscle soreness during vibration training is mostly due to accumulated fatigue over time or long term stress causing the body to tighten up. This condition can be improved through continuous whole body vibration training.

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PAIN 痛 KESAKITAN

During vibration training, one may feel pain on certain body parts. This is due to long-term accumulation of fatigue. This condition is more noticeable for those who do not exercise regularly, but continuous vertical vibration training will gradually enhance the body's circulation, and this condition will eventually go away.

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Whole Body Vertical Vibration Training (Passive Rope Skipping Exercise)

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Whole Body Vertical Vibration training is highly-efficient and covers areas from hairs to toes; nothing stays still, which is unattainable by ordinary exercises. Whole Body Vertical Vibration training produces a harmonious resonance frequency at around 5-13Hz, similar to the ones produced by QiGong, WaiTanGong, and ChanDouGong (ancient Chinese exercises). A research from the National Science Council with the title "There is Indeed Scientific Basis for WaiTanGong" has finally found the answer by looking for the "flow of Qi" of WaiTanGong through digital imaging technology. The research concluded that blood flow can be enhanced through vibration waves produced by WaiTanGong, and illnesses will be reduced naturally once the blood flow is smooth. The scientific principle of WaiTanGong's capability to eliminate illnesses and improve health is that the whole body will vibrate and show a resonance phenomenon after the WaiTanGong produces the "Qi", which is when the body vibrates at its own natural frequency and it can be as high as 5.6Hz (5.6 times per second). The body will eliminate illnesses and become healthy naturally because when in a relaxed state, blood flow resistance is reduced, blood flow is smooth, the cells in each body part are activated, each organ functions well, heart load is reduced, and heart rate is well regulated.

The notion that WaiTanGong is commonly known as an exercise of health and happiness has its scientific and medical justifications. The vibration exercise of WaiTanGong has unique effects on blood circulation, which other exercises cannot rival.

Every microscopical part of the body will be vibrating under the effect of the vibration. This eases the dilation of blood vessels and blood flow will accelerate naturally, which in turn improves microcirculation. This is entirely different from the ordinary exercise that causes muscle stiffness and impulsive impact. Normal organ function is dependent on how smooth the blood supply is. The human body, from hair to toe and from viscera to the skin, is formed by various kinds of cells, which is dependent on blood transportation for required nourishment and oxygen. The flow rate and volume of blood directly affect one's health; the quicker the blood flow and the higher the blood volume, the more oxygen and nutrient supply to the cells. As a result, the cells will be more active and organs will be healthier. One may even develop rosy skin complexion and become dexterous, with a healthy appearance. Hence, promoting blood circulation to supply sufficient oxygen and nutrients to the whole body to ensure normal organ function is the basic requirement and principle in health improvement.

From another point of view, every inch of the body muscle is vibrating during vibration, even the facial muscle. The spine will vibrate front and back in a wave-like motion. Every section of the body is involved in the lateral exercise. The front and back vibration of the entire spine exerts a massaging effect on the viscera, and the body gradually heats up and perspires. This kind of perspiration is due to the effect of vibration of each body part, which is different from the perspiration caused by muscle stiffness and forceful exertion. Therefore, you will feel a sense of well-being, your pores open up, and your body will be at ease after the preparatory vibration training. The mystery of vibration is insightful. From the research "There is Indeed Scientific Basis for WaiTanGong", it shows similarities between the theory of WaiTanGong and Whole Body Vertical Vibration training.

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Whole Body Vertical Vibration Training VS Exercise Regime for Visceral Improvement

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- Reprinted from Kai-Qiang Huang, Visceral Revolution -

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- Dicitak semula daripada Kai-Qiang Huang, Revolusi Visera -

“It Is Better to Go for a Walk than Training, It Is Better to Vibrate than Walking”, “vibration” is one of the traditional training methods of the Chinese. Due to the fact that modern people are experiencing greater pressure in their life and work, health care professionals recommend daily vertical vibration to eliminate body and mental fatigue, improve visceral vitality and to boost health. The body capillaries will be rejuvenated and blood circulation will be smoother due to body vibration. This can stave off issues like headache, poor digestive system, constipation, and muscle ache, and even improve immunity. ChanDouGong is not an intense exercise but a mild and comfortable exercise, which is not restricted by time, space and age.

Everyone knows that one need to move in order to be alive and exercise is an activity beneficial to health. But there are people who could not exercise persistently even if they are determined. There are also people who do various physical exercises but could not get any curative effect. All these are due to human idleness, which hinders people from getting sufficient appropriate exercise.

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Exercise regime for visceral improvement is simple and easy, and it is a health care method that can prevent and treat illnesses. However, the guiding principle is to not over exercise, overstrain the muscles and overexert you. The best exercise regime is to appropriately vibrate each organ so that the meridians can be smooth, and blood and "Qi" can circulate throughout the body.

Most people think that one needs to pant and perspire profusely in order to experience better results. The fact that the body will be exhausted due to overexertion if the exercise is too intense, which is counterproductive. The body will need a lot of rest and food to compensate, making it vulnerable to injuries due to ageing factor. The visceral arrangement is actually acclimatised for exercise. Humans are acclimatised for upright walking, running and jumping and that is why viscera are arranged upright, making it possible for energy to reach the viscera.

Visceral improving vertical vibration training regime emphasises on aerobic, mild and persistent exercise, which is to exercise for 30 minutes continuously each time, while maintaining a rosy complexion, perspiring but without being out of breath. Complementing the exercise with stretching exercises can give the limbs and spine a good stretch to open up the body's circulation. This will enhance metabolism and eliminate toxins from the cells, and achieve the effect of visceral recuperation.

Visceral improving exercise method is similar to Whole Body Vertical Vibration training promoted in many foreign countries. Both achieve body building and health care effects through vibration from the bottom to top.

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There are no side effects. Vibration is the safest health-preserving exercise. Due to individual differences, some may experience numbness, drowsiness, dizziness, fatigue, itchiness, aches, shoulder stiffness, increased frequency in passing gas and increased bowel movement after using the vibration device. These symptoms are normal, and will improve after a week of continuous use. If the symptom(s) do not go away after a week, then there may be a certain underlying illness and a health examination may be required.

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Q2 WHAT IS THE RECOMMENDED DURATION OF USE IN A DAY?

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Usage time is dependent on the user's goal and body condition. For general health, use once in the morning and once in the evening, 15 minutes each time. To strengthen the body, use once in the morning and once in the evening, 30 minutes each time. To tone up, use 3-4 times daily, 30 minutes each time. Usage may be repeated after 1 hour of rest after each session, 15-30 minutes each time. Frequency may be increased gradually once the body has adapted, and to obtain the best results, the frequency needs to be within the body's acceptable limits.

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Yes, it can. There are many successful cases where stroke patients experience speedy recovery in physical functions by utilising vertical vibration. In other words, vertical vibration is great for stroke patients. They may hold on to the edge of a bed, furniture and etc. for support. Frequent use will exert more health promoting benefits.

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Q4 HOW SHOULD THE DEVICE BE MAINTAINED?

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Limit it to indoor use and avoid places with high moisture. Please handle with care when moving the device, avoid dropping or collision to prevent damage to the device. Do not wipe the vibration device with volatile oil, alcohols and etc. as it might damage the device's surface. No special maintenance is required.

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Please start from a low speed when using the vibration machine; speed may be increased once the body has warmed up or once acclimatised. It works best on user's most comfortable speed level.

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HEALTH PRESERVING VIBRATION 基礎養生律動

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Vibration health care is so easy
律動保健好輕鬆

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There is no easier exercise than vibration; one would only need to stand or sit to eliminate illnesses and regain health without perspiring. Suitable for the whole family.

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10-15
 reps each time
 每次10-15循環
 10-15 ulangan
 setiap kali

Steps: Stand with feet shoulder-width apart, relax the body and take a deep breath and at the same time, open up palms and stretch arms backwards, and exhale slowly. While exhaling, slowly raise your arms, then slowly lower your arms.

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Steps: Stand with feet shoulder-width apart; bend elbows upward to form a 90° angle. Stretch out the upper arms to the sides while taking a deep breath and engage abs. Exhale and slowly raise up the arms at the same time. Lastly, drop and hold arms straight out to the sides.

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Standing Shoulder Squeeze
雙手後扣
Regangan Bahu Berdiri

10-15
reps each time
每次10-15循環
10-15 ulangan
setiap kali

Steps: Stand with feet shoulder-width apart. Relax the body and take a deep breath and at the same time, clasp hands behind the lower back. Gradually stretch the hands backward, relax and exhale slowly.

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10-15
reps each time
每次10-15循環
10-15 ulangan
setiap kali

Steps: Stand with feet shoulder-width apart. Clasp hands in front of the chest, inhale slowly and thrust the chest forward and contract the lower abdomen. Pull hands in opposite directions and exhale slowly at the same time.

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Steps: Kneel down on the floor; place hands shoulder-width apart on the device for support. Relax the neck and shoulders, and lower down the head with arms straight.

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Steps: Stand on the device with feet shoulder-width apart. Bend upper body forward, hold on to the back of a chair with both hands and relax with legs straightened.

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Steps: Stand on the device with feet shoulder-width apart. Slightly bend the upper body, place hands in front of the chest or hold on to wall, table or chair for support; relax and inhale slowly. Exhale slowly while contracting the anus and hips. Repeat.

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Leg Soothing 舒緩腿部 Melegakan Kaki

Steps: Sit on a chair and relax, place feet on the device. Alternatively, you may sit or lie on the floor and place calves on the device.

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Planking 平板支撐 Planking

Adjust the duration based on individual capacity

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Steps: Place elbows flat on the device; upper arms and platform will be at a 90° angle. Support the body with tiptoes and form a straight line from head to heels. Contract the lower abdomen and hips. Switch position of the arms and legs and repeat. Adjust planking time based on individual capacity.

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Pressing QingMing Acupoint
指按睛明穴
Menekan Titik Akupunktur QingMing

15
minutes each time
每次15分鐘
15 minit
setiap kali

Steps: Hold forehead with hands; gently press acupoint tissue around the eyes with thumbs.

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15
minutes each time
每次15分鐘
15 minit
setiap kali

Steps: Sit on a chair and place feet on the device; hold chin with both hands.

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Acupoint Massage
按摩穴位
Urutan Titik Akupunktur

15
minutes each time
每次15分鐘
15 minit
setiap kali

Steps: Sit on a chair and place feet on the device, gently massage temples with thumbs.

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Stretching is one of the health preserving methods unanimously recommended by health care professionals globally. Stretching does not only limber the body, but it also promotes and preserves health. It is highly beneficial to health for muscles and bones remain nimble.

As the saying goes "it is difficult to have a long life with rigid old muscles", and from it we learn how important stretching is to humans.

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Stretch 5-10 seconds per movement for 5 minutes

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Steps: Stand with feet shoulder-width apart; lift arms horizontally to the sides. While taking a deep breath, slowly lift arms upward. Exhale slowly and clasp palms.

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Steps: Stand with feet shoulder-width apart; lift arms horizontally to the front. With palm facing to the front, pull right fingers downward with left hand. Switch sides and repeat.

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Triceps Stretch 手臂後扣拉筋 Regangan Triceps

Stretch 10 seconds
per movement for
5 minutes

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Regangan 10 saat
bagi setiap pergerakan
selama 5 minit

Steps: Stand feet shoulder-width apart. With arms extended overhead, bend left elbow and reach your right hand overhead and grasp just below your left elbow. Gently pull your left elbow down and toward your head. Switch sides and repeat.

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- Stretches arm muscles, boosts body metabolism

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Regangan 10 saat
bagi setiap pergerakan
selama 5 minit

Steps: Stand with feet shoulder-width apart. Lift and straighten the right arm to the left, with lower body remaining still, bend the left elbow to hook the right arm. Switch sides and repeat.

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- Relieves soreness and enhances mobility of the shoulder

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2

3

Steps: Lift head slightly by extending the neck upwards to form a 15-20° angle. Relax the shoulders and stand on the device. Slowly roll the shoulders forward, up, back and down, in clockwise direction.

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2

3

Steps: Lift head slightly by extending the neck upwards to form a 15-20° angle. Relax the shoulders and stand with feet shoulder-width apart. Relax the body, take a deep breath slowly, clasp hands behind the back and extend shoulders backward gradually. Bend forward and extend clasped hands as far as possible, stretching the legs, and exhale slowly.

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Steps: Kneel down on all fours on the device. Gently lower the back towards the ground while inhaling with limbs remaining still. Arch the back upwards as much as possible, relax and exhale slowly.

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Arm Cross Stretching Behind the Back
雙手後勾交扣

Regangan Tangan Melintang Bahagian Belakang

Stretch 3-5 seconds per movement for 5 minutes

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Steps: Reach right arm straight up, then bend elbow and let hand fall behind your head. Move left arm behind the back and bend the arm, letting the back of your left hand rest against the right shoulder blade. Reach to grab right fingertips with the left hand pulling into opposite directions. Repeat on the other side.

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Steps: Kneel down on all fours on the device. With the lower body remaining still, stretch the side of waist gently to the left leading with the head. Switch sides and repeat.

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Steps: Stand with feet shoulder-width apart, clasp hands in front of the chest, with the lower body remaining still. Take a deep breath, keep the abdominal muscles tight and rotate the upper body to the right in 2 stages.

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* It is important to remember that lower body has to remain still to stretch both sides of the side waist muscles to its full range

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Steps: Stand with feet shoulder-width apart; lift the arms horizontally to the side. Slowly bend the waist to the left; stretch as much as possible with right arm over the head. Switch sides and repeat.

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Steps: Put left foot on the platform and bend knees at around 90°; place hands on thigh, slightly tilt body to the front, keep the abdominal muscles tight, straighten and extend your right leg out to the right. Switch sides and repeat.

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Steps: Stand next to the device, slightly bend the right leg; straighten the left leg to the front and hook the heel on the platform. Bend the body forward, straighten the back, and hold onto thighs. Switch sides and repeat.

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Steps: Sit on the device, with the soles of your feet together, knees bent out to sides. Hold onto your feet, engage abs, and slowly lower body towards feet as far as possible while pressing knees toward the floor. Hold for 3-5 seconds before relaxing. Repeat this movement continuously.

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Steps: Put right foot onto the platform; bend right knee to form a 90° angle. Extend left leg behind, slightly bend and tiptoe on the floor. With the upper body remaining straight, engage abs. Lean upper body forward, place hands on right knee, push and straighten to extend the left leg backward. Switch sides and repeat.

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Steps: Sit on the device; extend right leg out to the right, with toes facing upward. Tuck left foot to inner right thigh. Tilt the head forward. Extend right hand as far as possible towards the toes of the right foot and hold it. Switch sides and repeat.

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Steps: Stand with feet shoulder-width apart; lift the arms horizontally to the front. Squat while exhaling slowly. While engaging abdominal muscles, lower squat position with hands and knees remain still. The joint angle of the inner legs should not be lesser than 90°.

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Steps: Sit on the device, cross your left leg over your right, and place your left foot flat on the platform. Rotate upper body to the left while the lower body remains still, keeping the back straightened. Switch sides and repeat.

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MUSCLE STRENGTH AND BODY SCULPTING TRAINING

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The secret of having strong muscle that you must know

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The stronger the muscle, the better the microcirculation

The stronger the muscle, the faster the reflux rate of venous blood

The stronger the muscle, the lesser the load on joints

The stronger the muscle, the higher the body temperature, the stronger the immunity

The stronger the muscle, the younger the organs

The stronger the muscle, the lower the body fat

The stronger the muscle, the better the figure and confidence, the better the spirit

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Steps: Stand with feet shoulder-width apart. Hold and lift dumbbells horizontally to the front. Pull arms back to the sides to expand the chest while exhaling. Take a deep breath, lift and cross the arms above the head. While exhaling, lower arms to form an 'O' shape and then a 'U' shape.

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Steps: Kneel on the floor and cross calves; straighten arms and place on the device, engage abs, bend elbows and lower down upper body while exhaling.

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Steps: Form a push-up position with arms straight on the device. Ground toes into the floor, position body in a straight line, bend the elbows to lower the body while exhaling.

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Planking
棒式支撐
Planking

Steps: Ground toes on the device as support and position body in a straight line; bend elbows to 90° to support the body, lower the body while engaging core muscles. Switch positions between the arms and legs and repeat.

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Steps: Sit on the device with knees bent and feet firmly planted on the floor, bend right elbow placing fingers behind right ear; rotate waist to the right and straighten left arm to touch right knee. Switch sides and repeat.

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✦ Put a towel, pad or yoga mat underneath the elbows to alleviate pressure when placing elbows on the device
手肘置於機台時，可鋪放毛巾、軟墊或瑜珈墊來減緩手肘的不適壓力

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Regangan 3-5 saat bagi setiap pergerakan selama 5 minit

Adjust posture based on individual capacity

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Steps: Sit on the device with knees bent and feet firmly planted on the floor; reach back and hold onto the sides of the device, lift right leg so that the calf is parallel to the floor while inhaling; exhale slowly, straighten and kick right leg out while engaging the abs. Switch sides and repeat.

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Adjust posture
 based on individual
 capacity

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Laraskan postur getaran
 berdasarkan kapasiti
 individu

Steps: Sit on the device with knees bent and feet firmly planted on the floor, lean back and form a 30° angle with the device, while holding arms out to the front, engage abs while exhaling; use the strength of the abs to move the upper body forward to the original sitting position with feet firmly planted on the floor.

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Sit-up with Body Rotation
仰臥轉體
Bangkit Tubi dengan Pusingan Badan

Advanced Movement
進階動作
Senaman Tahap Sukar

Adjust posture based on individual capacity
依自身能力調整律動姿勢

Laraskan postur getaran berdasarkan kapasiti individu

Stretch 3-5 seconds per movement for 5 minutes
每一動作伸展約3-5秒持續5分鐘

Regangan 3-5 saat bagi setiap pergerakan selama 5 minit

Steps: Sit on the device, bend and lift knees parallel to chest level, place hands supporting the head and engage abs; rotate upper body to the left while exhaling; move left leg closer to the chest with left knee touching the right elbow, switch sides and rotate the upper body to the right; straighten and extend the right leg out, switch sides again - completing 1 rep.

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Laraskan postur getaran berdasarkan kapasiti individu

Put a towel, pad or yoga mat underneath the elbows to alleviate pressure when placing elbows on the device
 手肘置於機台時，可鋪放毛巾、軟墊或瑜珈墊來減緩手肘的不適壓力

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Steps: Position the body in a straight line by supporting the body with right hand extended on the device; put feet together, straighten and lift left arm up high parallel with the right hand; keep abs engaged, hips remain still, then lift and stretch the left leg. Switch sides and repeat.

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Steps: Support the body with right forearm on the device and position the body in a straight line; left arm stays on left side, put feet together and engage abs; use side waist strength to lift hip upward with elbow remaining still; lift the left arm up high. Switch sides and repeat.

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Plank Roll-in
瑜珈球捲腹
Plank Mengguling

Steps: A. Place yoga ball and sit on it; lift arms horizontally as a relaxing movement. B. Support body by placing straightened hands on the device; place calves on the yoga ball, while engaging abs; tuck legs and roll the yoga ball towards the abs, then slowly back to Position 1.

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Steps: Sit on the device with sit-up position, lift the legs to a 30° angle off the floor, bend elbows and hold on to the sides of the device; place yoga ball on the concave abdomen, engage abs and slowly lower down legs while exhaling without thighs touching the platform.

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Deep Squat in Stages 分段深蹲 Cangkung Rendah-rendah Secara Berperingkat

Stretch 5-10 seconds
per movement for
5 minutes

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Regangan 5-10 saat
bagi setiap pergerakan
selama 5 minit

Steps: Stand with feet shoulder-width apart, lift arms horizontally to the front and lower squat in stages while exhaling; engage abs, with arms and knees remaining still, squat down slowly until thighs are parallel with the device. Bend elbows, while trying your best to tiptoe, remaining stable and balanced.

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Headache, Sore Eyes
頭痛、眼睛酸澀
Sakit Kepala, Kelelahan Mata

Vibrate 30 seconds
for each movement
Repeat Training

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+

Vibrate 30 seconds for each movement for 15 minutes, at least twice a day, once in morning and once in the evening. Vibration can be repeated after every 2 hours of rest for better effect.

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2

3

Vibrate 15 minutes each time, at least once in morning and once in the evening every day. Repeat training after every 2 hours of rest for better effect.

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Vibrate 15-30 minutes each time, at least once in morning and once in the evening every day. Repeat training after every 2 hours of rest for better effect.

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* For those who are suffering from both osteoporosis and degenerative arthritis, are recommended to start by sitting on the chair and with feet on the vibration device, knees bent at an angle of more than 90° to improve the arthritis pain effectively. When the knee joints are no longer painful, use it in standing posture or adopt the posture as in Picture B, which relaxes the lower back muscles and increases bone density at the same time.

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
Getaran untuk Pemulihan Badan

Health Benefits of Vibration
律動對於健康的好處
Manfaat-manfaat Kesehatan Getaran

Vertical Vibration has tremendous benefits in preventive health care, health recuperation and anti-ageing. Decades of research literatures have proven that other than improving bone density, preventing osteoporosis, increasing muscle strength and balancing, Vertical Vibration also assists in recovery of bodily functions or improvement of the degenerative issues, ageing problems, modern sub-health problems and etc.

Vertical Vibration has gone beyond sport and medical level and become a non-medical and non-invasive supplementary health promoting method in the medical field after decades of evolution. Vertical Vibration is the health promoting exercise for the whole family, anytime and anywhere.

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Visceral Prolapse Prevention
內臟下垂
Pencegahan Prolaps Visera

Adjust the duration
based on individual
capacity

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Vibrate 5 minutes continuously each time, at least once in the morning and once in the evening every day. Repeat after every 2 hours of rest for better effect.

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Vibrate 3-5 seconds per movement for 5 minutes, at least once in morning and once in the evening every day. Repeat training after every 2 hours of rest for better effect.

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Research of Vertical Vibration Therapy by Experts and Scholars

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Vibration Suppresses Production of Body Fat

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Getaran Membantu Penghasilan Lemak Badan

Conducted research in reducing body fat by using Whole Body Vertical Vibration device, discovered that Whole Body Vertical Vibration can inhibit obesity caused by high fat diet, prevent health problems associated to excess body fat accumulation or obesity. It is an effective health-promoting and preventive measure.

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Vibration Promotes Energy Consumption of the Body

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Getaran Mendorong Penggunaan Tenaga Badan

Conducted research in energy consumption by using Whole Body Vertical Vibration device, results indicates that Vertical Vibration significantly increases energy consumption and is beneficial for overweight people.

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Conducted research by using Whole Body Vertical Vibration device, results indicates that Whole Body Vertical Vibration boosts the vitalisation of muscle cells and quadruples the collagen multiplication to firm up skin.

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Jabatan Sains Sukan, Universiti Sukan Nasional Taiwan

Whole-Body Vibration Training Effect on Physical Performance and Obesity in Mice

Chi-Chang Huang, Tzu-Ling Tseng, Wen-Ching Huang, Yi-Hsiu Chung,

Hsiao-Li Chuang, and Jyh-Hong Wu

International Journal of Medical Sciences

2014; 11(12): 1218-1227. doi: 10.7150/ijms.9975

Professor Soun-Cheng Wang

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Professor Soun-Cheng Wang

Professor in Graduate Institution of Sport and Leisure Education,

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Professor di Institut Pendidikan Sukan dan Riadah Siswazah,

Universiti Nasional Chung Cheng

The Influence of Whole Body Vibration Stimulus on Energy Expenditure in Obesities

Wang, Soun Cheng Ph.D.; Chiang, Tai Ju; July, 2007

Master's Thesis Graduate Institute of Sport and Leisure Education ;

National Chung Cheng University

Medical Professor Chia-Hsin Chen

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Profesor Perubatan Chia-Hsin Chen

Director of Rehabilitation Center, Kaohsiung Medical University

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Low-magnitude vertical vibration promoting myotube formation in myoblasts

Chia-Hsin Chen, M.D Department of Physical and Rehabilitation Medicine,

Kaohsiung Medical University Hospital, College of Medicine

Promotes Proliferation and Repair of Nerve Cells

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Japan conducted research by using Whole Body Vertical Vibration device and discovered that low-frequency Whole Body Vertical Vibration can help in the proliferation of nerve cell and recovery of damaged nerve cell.

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Increases Bone Density

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Meningkatkan Ketumpatan Tulang

Research conducted on bone density improvement by using Whole Body Vertical Vibration device for 5 days per week, 30 minutes every day; showing significant increase in bone density after 5 weeks of training.

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Induces Recovery from Fatigue and Lactic Acid Metabolism

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Conducted research on fatigue recovery after exercise by using Whole Body Vertical Vibration device, results indicates that Whole Body Vertical Vibration has significant effects on fatigue recovery after exercise and lactic acid metabolism.

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Conducted elderly vibration training for 6 weeks by using Whole Body Vertical Vibration, showing significant benefits on arteriosclerosis.

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Professor Tutida R

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Professor Tutida R

Department of Occupational Therapy cum Faculty of Health and Welfare, Hiroshima University

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Low-Frequency, Whole Body Vibration Induced Neurite Outgrowth by P12m3 Cells with Impaired Nerve Growth Factor-Induced Neurite Outgrowth
Yoshihisa Koike, Reiko Tutida, Yuko Hayashi, Yoko Yamanishi and Yoshio Kano
J Nov Physiother ISSN:2165-7025 JNP, an open access journal Volume 5 • Issue 1 • 1000249

Medical Professor Chia-Hsin Chen

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Professor Perubatan Chia-Hsin Chen

Director of Rehabilitation Center, Kaohsiung Medical University

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Whole Body Vibration Increasing Bone Density in the Osteoporotic Rats by Tail Suspension
Chia-Hsin Chen, Gwo-Jaw Wang, Mei-Ling Ho, Mao-Hsiung Huang, Chin Hsu.

Professor Ching-Feng Cheng

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Professor Ching-Feng Cheng

National Taiwan Normal University Department

of Athletic Performance Professor

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Professor Jabatan Prestasi Keolahragaan, Universiti Normal Nasional Taiwan

Acute Effect of Vibration Exercise on Grip Strength and Physiological Recovery from High-Intensity Intermittent Grip Exercise
Ching-Feng Cheng; Chia-Lun Lee; Yi-Shan Yang (2010/09)

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Professor Sains Sukan, Profesor Chen, W.C.

Director of Sports Room, Chang Gung University of Science and Technology

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Influence of Six Weeks Whole-Body Vibration Exercise On Arterial Stiffness In The Elderly 16th Annual Congress Of The European College Of Sport Science (Ecsc)
Chen, K.W.C.1, Chen, W.C., Chia, P.S., Chi, S.F., Ye, J.Y., Wang, S.W.

Increases the Concentration of Nitric Oxide

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Meningkatkan Kepekatan Nitrik Oksida

"After 6 weeks of Whole Body Vertical Vibration training, systolic blood pressure had an average decrease of 16mmHg, diastolic blood pressure had an average decrease of 6.4mmHg, pulse pressure had decreased 12.2mmHg, and nitric oxide (NO) concentration had greatly increased in the body."

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Helps Preserve Body Tissue to Its Youthful Best

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Membantu Meremajakan Tisu-tisu Badan Sebaik Mungkin

Conducted study using Whole Body Vertical Vibration device and found that vibration can enhance the functionality of bones, muscles, liver and kidney tissues, increase sports performance, eliminate fatigue, and promote anti-ageing effects.

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Stimulates and Enhances Explosive Power and Running Efficiency of Runners

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Conducted research on scientific exercise training for track and field athletes with Whole Body Vertical Vibration device, and results showed a significant increase of explosive power of the runners' lower limbs and their running efficacy after short term high intensity vibration training.

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Influence of Six Weeks Whole-Body Vibration Exercise On Arterial Stiffness In The Elderly
Chen, K.W.C., Chen, W.C., Chia, P.S., Chi, S.F., Ye, J.Y., Wang, S.W.

16th Annual Congress Of The European College Of Sport Science (Ecss) (2011)

Professor Chi-Chang Huang

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Profesor Chi-Chang Huang

National Taiwan Sport University Department of Sports Science

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Jabatan Sains Sukan, Universiti Sukan Nasional Taiwan

Effect of whole-body vibration training on body composition, exercise performance and biochemical responses in aged mice

Corresponding Author: Prof. Chi-Chang Huang, Ph.D.

Lin CI, Huang WC, Chen WC, Kan NW, Wei L, Chiu YS, Huang CC.

Metabolism. 2015 Sep;64(9):1146-56. doi: 10.1016/j.metabol.2015.05.007

Professor Ching-Feng Cheng

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Profesor Ching-Feng Cheng

National Taiwan Normal University Department

of Athletic Performance Professor

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Profesor Jabatan Prestasi Keolahragaan, Universiti Normal Nasional Taiwan

Improvement in Running Economy After 8 Weeks of Whole-body Vibration Training
Ching-Feng Cheng, Kuo-Hui Cheng, Yu-Ming Lee, Hsin-Wei Huang,

Yu-Hsuan Kuo, Heng-Ju Lee

Journal of Strength and Conditioning Research Publish Ahead of Print DOI:

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Conducted vibration therapy research on stroke patient by using Whole Body Vertical Vibration device, which indicated Whole Body Vibration had significant improvement effects on cramps and mobility of stroke patient.

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Increases Functional Fitness in the Elderly

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Conducted research on elderly mobility by using Whole Body Vertical Vibration, the results indicates that Whole Body Vertical Vibration is a safe and scientifically-proven exercise method that significantly boosts elderly's mobility.

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Relieves Constipation

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Melegakan Masalah Sembelit

Conducted research on constipation improvement by using Whole Body Vertical Vibration device, results indicates that Whole Body Vertical Vibration has significant effects on improving constipation.

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Enhances Sports Performance and Fights Fatigue

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Conducted a 4-week study on scientific training with Whole Body Vertical Vibration device and confirmed its benefits on enhancing sports performance and fighting fatigue, without causing physiological and biochemical side effects. The use of Whole Body Vibration training in sports training has both health promoting and sport safety properties.

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Effects of a single session of whole body vibration on ankle plantar flexion spasticity and gait performance in patients with chronic stroke: a randomized controlled trial
Chan KS, Liu CW, Chen TW, Weng MC, Huang MH, Chen CH.
Clin Rehabil. 2012 Dec;26(12):1087-95. doi: 10.1177/0269215512446314

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Profesor Sains Sukan, Profesor Chen, W.C.

Director of Sports Room, Chang Gung University

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Effect of 6 weeks whole-body vibration exercise on functional fitness in the elderly
CHEN, K. W. C., CHEN, W. C., CHIA, P. S., CHIU, P. K., WANG, S. W.

Doctor Tsung-Ju Wu

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Doktor Tsung-Ju Wu

Director of Department of Rehabilitation, Caotun Branch,

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Pengarah Jabatan Pemulihan, Cawangan Caotun, Hospital Kristian Changhua

The effect of whole-body vibration on improving the symptoms of functional constipation: A randomized controlled pilot study
Tsung-Ju Wu, Ta-Sen Wei, Ya-Hui Chou, Chao-Pin Yang, Chien-Lin Wu, en-Chun Chen, and Sen-Yung Liu

Professor Chi-Chang Huang

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Profesor Chi-Chang Huang

National Taiwan Sport University Department of Sports Science

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Jabatan Sains Sukan, Universiti Sukan Nasional Taiwan

Whole-Body Vibration Training Effect on Physical Performance and Obesity in Mice
Chi-Chang Huang, Tzu-Ling Tseng, Wen-Ching Huang, Yi-Hsiu Chung, Hsiao-Li Chuang, and Jyh-Hong Wu
International Journal of Medical Sciences
2014; 11(12): 1218-1227. doi: 10.7150/ijms.9975

Health Exercise Modes

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Mod Latihan Kesihatan

7 Fixed Frequency Selection

7 段固定頻率調整
7 Pelarasan Frekuensi Tetap

- 1 General Health Vibration**
舒活養生律動
Getaran Penjagaan Kesihatan Am
- 2 Leg Muscle Relaxation Vibration**
腿肌舒壓律動
Getaran Melegakan Otot Kaki
- 3 Waist and Abdomen
Massaging Vibration**
腰腹按摩律動
Getaran Mengurut Pinggang & Abdomen
- 4 Leg and Hip Vibrating Vibration**
腿臀彈抖律動
Getaran Kaki dan Pinggul
- 5 Shoulder and Neck Vitalisation
and Relaxation**
肩頸舒緩放鬆
Pencergasan & Pelegaan Bahu dan Leher
- 6 Whole Body Muscle Vitalisation**
全身肌肉舒活
Pencergasan Otot Seluruh Badan
- 7 Whole Body Deep Stimulation**
全身深層刺激
Rangsangan Mendalam Seluruh Badan

3 Preset Automatic Mode Functions

3 組自動程序功能
3 Fungsi Mod Automatik Pra-tetap

- M1 Alleviation and Relaxation
Recuperation**
舒緩放鬆調理
Pemulihan Bersifat Relaksasi dan Melegakan
- M2 Massage and Vitalisation**
按摩活化調理
Pemulihan Mengurut dan Mencergaskan
- M3 Internal Circulation Recovery**
體內循環調理
Pemulihan Peredaran Dalam